

SP8-U-FLIM

Pos.	Objective	DIC	Magnification	Num. Aperture	Imm.-Medium	Working distance (mm)	Internal number
1	HCPL FLUOTAR 10x/0,30 (506505)	-	10	0,30	air	11,0	245
2	HC PL APO CS2 IMM CORR 20x/0,75 (506343)	+	20	0,75	water, oil, gly	0,67	246
3	HC FLUOTAR L VIS IR 25x/0,95 W (506374)	+	25	0,95	water	2,50	247
4	HC PL APO CS2 40x/1,10 W (506357)	+	40	1,10	water	0,65	248
5	HC PL APO CS2 63x/1,20 W CORR (506346)	+	63	1,20	water	0,30	249
6	HC PL APO CS2 63x/1,40 OIL	+	63	1,40	oil	0,14	250
7	HC PL APO CS2 40x/1,30 OIL	-	40	1,30	oil	0,24	256

SP8-U-FLIM

Filters for direct visualization by eye

Filter Pos.	Filterset	Fluorochrome	Filter number
1	empty		-
2	38 HE GFP AL-	GFP	
3	43 HE AL-	DsRED, RFP	
4	49 DAPI AL-	DAPI	
5	C-DIC RL AL-	DIC	
6			
for exchange:			